

EMG-triggered surface FES for arm reaching in tetraplegia

Dimitra Blana, Edward K Chadwick

Institute for Science and Technology in Medicine, Keele University



Neil Postans, Simon Pickard



Robert Jones and Agnes Hunt Orthopaedic Hospital, Oswestry




Implanted FES for restoring arm function



Implanted FES for restoring arm function




Surface FES for restoring arm function?



Surface FES




Orthosis



Model-based
FES controller




Surface FES for restoring arm function?



Surface FES



Orthosis




**Model-based
FES controller**



How do we use models?

- To optimise the design



How do we use models?

- To optimise the design
- To be a testing platform



How do we use models?

- To optimise the design
- To be a testing platform
- To act as (part of) the device controller




How do we use models?

- To optimise the design
- **To be a testing platform**
- To act as (part of) the device controller



EMG-triggered surface FES

The patient was a 56-year old male, 19 months post injury (C4 Frankel A with some denervation at C8, but stimlatable C5,6 & 7).



EMG-triggered surface FES

The patient was a 56-year old male, 19 months post injury (C4 Frankel A with some denervation at C8, but stimlatable C5,6 & 7).

We used an Ottobock STIWELL surface stimulator and placed recording EMG electrodes on the biceps and anterior deltoid, and stimulating electrodes on the triceps and wrist extensors.




The algorithm

Reaching out:



Bringing the arm in:




Example computer simulation

The deltoid activity is slowly increased while the biceps is off.

When the deltoid reaches the threshold of 10%

the triceps is stimulated and the elbow fully extends.




Example computer simulation

The deltoid activity is slowly increased while the biceps is off.

When the deltoid reaches the threshold of 10%

the triceps is stimulated and the elbow fully extends.




Example computer simulation

The deltoid activity is slowly increased while the biceps is off.

When the deltoid reaches the threshold of 10%


the triceps is stimulated and the elbow fully extends.



Example computer simulation

The biceps activity is then slowly increased.


When it reaches a threshold of 5% the triceps stimulation is turned off and the elbow flexes.



Example computer simulation

The biceps activity is then slowly increased.


When it reaches a threshold of 5% the triceps stimulation is turned off and the elbow flexes.



Graphical Interface

Sliders for deltoid
and biceps “EMG”

control the
triceps stimulation



Next steps

EMG-controlled

- Test ~~EMG-triggered~~ surface FES for arm reaching in tetraplegia
- Explore other technologies that could work with FES
- Optimize the choice and combination of EMG control signals using personalised computer simulations

